

# Development and characterization of 3D semiconductor X-rays detectors for medical imaging

Marie-Laure Avenel, Eric Gros d'Aillon

CEA-LETI, DETectors Laboratory

marie-laure.avenel@cea.fr

# Outlines


- Problematic
- 3D detector process
- 3D CdTe detector
- First results for 3D GaAs detector

## Study context and problematic

- **Application:** Radiography X-rays detectors (20-60 keV energy range) using
  - A semiconductor material
  - A micro-structured detection geometry

# Study context and problematic

- **Application:** Radiography X-rays detectors (20-60 keV energy range) using
  - A semiconductor material
  - A micro-structured detection geometry
- **Problematic :**



## 3D detector process



Bulk Semiconductor


# 3D detector process



Bulk Semiconductor





Laser drilling\* by a Nd-YLF laser  
(@263 nm, 40 ns, 600 μJ, 3 kHz)  
in percussion mode



\*D. Farcage, CEA-Saclay, LILM

# 3D detector process



\*D. Farcage, CEA-Saclay, LILM


# 3D detector process



Bulk Semiconductor



Laser drilling\* by a Nd-YLF laser  
(@263 nm, 40 ns, 600  $\mu$ J, 3 kHz)  
in percussion mode





\*D. Farcage, CEA-Saclay, LILM



Chemical etch to clean  
surface and holes




Aspect ratio until 50 : 1




Metal electrodes  
deposition by electroless  
or sputtering



# 3D detector process




# 3D detector process





# 3D CdTe detector

## Sample



- CdTe (Cl) Acrorad ; thickness = 1,6 mm
- 3x3 holes matrix ( $\Phi$  entrance = 120 /  $\Phi$  exit = 40  $\mu\text{m}$  ; pitch = 350  $\mu\text{m}$ )
- Electrodes : gold electroless + gold wires ( $\Phi$  = 25  $\mu\text{m}$ )



## Biasing configurations :




## Measurement under $^{241}\text{Am}$ and $^{57}\text{Co}$ gamma ray irradiation



## Measurement under synchrotron X-ray irradiation

- Mapping the detection response : measurement of charge efficiency
- Comparison with charge efficiency simulation

# 3D CdTe detector



- Demonstration of photon-counting ability using a non-optimized design (electrode diameter and pitch have to be optimized)
- Energy discrimination for both biasing configurations
- FWHM = 7 keV @ 60 keV and 16 keV @ 122 keV
- Spectrum tails due to inhomogeneous regions of charge collection and fluorescence in CdTe

## 3D CdTe detector

Synchrotron beam (ID6 beamline at ESRF) :

- Energy : 60 keV (radiography energy)
- Size : 20 x 20  $\mu\text{m}^2$
- Input rate :  $10^5$  photons/s

Scan parameters :

- 32 x 32-position top-surface scan of the detector at a pitch of 25  $\mu\text{m}$
- Spectrum collected for 10 s



# 3D CdTe detector

Synchrotron beam (ID6 beamline at ESRF) :

- Energy : 60 keV (radiography energy)
- Size : 20 x 20  $\mu\text{m}^2$
- Input rate :  $10^5$  photons/s

Scan parameters :

- 32 x 32-position top-surface scan of the detector at a pitch of 25  $\mu\text{m}$
- Spectrum collected for 10 s




Incident photon absorption efficiency is homogeneous over the whole surface of the detector

# 3D CdTe detector

Measurement of the absorbed charge : normalized peak channel (calibrated measurement set-up)

Measurement (-20 V)



The normalized peak channel is greater than 0.8 → good collection for both carriers


**The detector response is homogeneous**



# 3D CdTe detector

## Measurement of the absorbed charge and comparison with simulated Charge Induction Efficiency

Measurement (-20 V)







# 3D CdTe detector

## Measurement of the absorbed charge and comparison with simulated Charge Induction Efficiency

Measurement (-20 V)



Simulation (-20 V)





M. Ruat et al., "3D Semiconductor radiation detectors for medical imaging: simulation and design", *IEEE Nucl. Sci. Symposium Conference Record*, 2009

# 3D CdTe detector

## Measurement of the absorbed charge and comparison with simulated Charge Induction Efficiency

Measurement (-20 V)

Simulation (-20 V)



Four-leaf clover shape around central anode in experimental and simulated scan due to electric field distribution in biased 3D structure



M. Ruat et al., "3D Semiconductor radiation detectors for medical imaging: simulation and design", *IEEE Nucl. Sci. Symposium Conference Record*, 2009

# 3D CdTe detector

## Measurement of the absorbed charge and comparison with simulated Charge Induction Efficiency

Measurement (-20 V)

Simulation (-20 V)



Highest values of charge efficiency around cathodes due to transport properties of electrons upper the holes transport properties in CdTe


M. Ruat et al., "3D Semiconductor radiation detectors for medical imaging: simulation and design", *IEEE Nucl. Sci. Symposium Conference Record*, 2009




# 3D CdTe detector

## Measurement of the absorbed charge and comparison with simulated Charge Induction Efficiency

Measurement (-20 V)



Simulation (-20 V)




Low electric field and diffusion-dominated regions in the corner of square-shaped configuration: enhanced trapping and reduced charge collection efficiency

M. Ruat et al., "3D Semiconductor radiation detectors for medical imaging: simulation and design", *IEEE Nucl. Sci. Symposium Conference Record*, 2009


# First results for 3D GaAs detector

## GaAs 3D prototype :

- Semi-insulating GaAs Freiberger ; thickness = 600  $\mu\text{m}$
- 3x3 holes matrix :
  - Entrance diameter = 100  $\mu\text{m}$  / Exit diameter = 50  $\mu\text{m}$
  - Pitch: 150  $\mu\text{m}$
- Electrodes : Gold electroless



Results under gamma rays irradiation ( $^{241}\text{Am}$  et  $^{57}\text{Co}$ ) : Square-shaped biasing configuration @ -60 V



Demonstration of photon-counting ability

Efficiency for 60 keV = 2%

# Conclusion

## **3D geometry :**

Alternative detector for X-rays medical imaging applications for semiconductor whose transport properties are lower (than CdTe for example)

Requires technical effort

## **3D CdTe detector :**

Proof of concept for 3D geometry

Photon-counting ability

60 keV and 122 keV energy discrimination

Homogeneous absorption and charge collection under X-rays irradiation

## **3D Semi-insulating GaAs detector :**

First prototype has been developed

Photons counting under gamma rays irradiation

Perspectives :

- Electrode : Ti/Pt/Au sputtering
- Chemical etch or Reactive Ion Etching to clean holes surface
- Holes filling with metal and bonding

# leti

LABORATOIRE D'ÉLECTRONIQUE  
ET DE TECHNOLOGIES  
DE L'INFORMATION

CEA-Leti  
MINATEC Campus, 17 rue des Martyrs  
38054 GRENOBLE Cedex 9  
Tel. +33 4 38 78 36 25

[www.leti.fr](http://www.leti.fr)



# Thank you for your attention




énergie atomique • énergies alternatives



# Contribution of 3D geometry for GaAs

	CdZnTe (single crystal)	GaAs (epitaxy*)
$(\mu\tau)_e$ (cm <sup>2</sup> /V)	$3 \cdot 10^{-3}$	$8 \cdot 10^{-5}$
$(\mu\tau)_h$ (cm <sup>2</sup> /V)	$5 \cdot 10^{-6}$	$4 \cdot 10^{-5}$

\*Epitaxy GaAs: G.C. Sun *et al.*, NIM-A, 546, 2005





Hypothesis: full depletion



# Contribution of 3D geometry for GaAs

	CdZnTe (single crystal)	GaAs (epitaxy*)
$(\mu\tau)_e$ (cm <sup>2</sup> /V)	$3 \cdot 10^{-3}$	$8 \cdot 10^{-5}$
$(\mu\tau)_h$ (cm <sup>2</sup> /V)	$5 \cdot 10^{-6}$	$4 \cdot 10^{-5}$

\*Epitaxy GaAs: G.C. Sun *et al.*, NIM-A, 546, 2005




Hypothesis: full depletion


# Contribution of 3D geometry for GaAs

	CdZnTe (single crystal)	GaAs (epitaxy*)
$(\mu\tau)_e$ (cm <sup>2</sup> /V)	$3 \cdot 10^{-3}$	$8 \cdot 10^{-5}$
$(\mu\tau)_h$ (cm <sup>2</sup> /V)	$5 \cdot 10^{-6}$	$4 \cdot 10^{-5}$

\*Epitaxy GaAs: G.C. Sun *et al.*, NIM-A, 546, 2005



- Thickness
- CdTe 640μm
  - GaAs 2.3mm
- Voltage :
- 2/7, 1V
  - 10/50, 5V
  - 70/350, 50V
  - Planar CdTe 200V
  - Planar GaAs 720V




Hypothesis: full depletion


# Contribution of 3D geometry for GaAs

	CdZnTe (single crystal)	GaAs (epitaxy*)
$(\mu\tau)_e$ (cm <sup>2</sup> /V)	$3.10^{-3}$	$8.10^{-5}$
$(\mu\tau)_h$ (cm <sup>2</sup> /V)	$5.10^{-6}$	$4.10^{-5}$

\*Epitaxy GaAs: G.C. Sun *et al.*, NIM-A, 546, 2005



- Thickness
- CdTe 640μm
  - GaAs 2.3mm
- Voltage :
- 2/7, 1V
  - 10/50, 5V
  - 70/350, 50V
  - Planar CdTe 200V
  - Planar GaAs 720V



Structuring detector can make spectrometric a material whose transport properties are lower than CdZnTe

# Characterization of laser drilling

## Available Techniques :

- Dry etching
  - Chemical etching
  - Mechanical drilling
  - **Laser drilling**
- } Mainly developed for Silicon

## Advantages / disadvantages laser drilling :

- + : High aspect ratio
- + : Material independent
- : Damaged holes sides
- : Holes tapering
- : Serial process


## Laser drilling system : (CEA-Saclay)

A Nd-YLF pulsed laser has been used, exhibiting the following parameters:

- Wavelength = 263 nm, enabling a small spot size and small hole diameter
- Pulse energy = 600  $\mu$ J
- Rate = 3 kHz
- Pulse width = 40 ns

Holes are realised by laser beam percussion and samples are placed in the air.

Percussion  
drilling





# Characterization of laser drilling

## Feasibility of laser drilling

- 60 000 holes matrix drilled into 1 mm of CdTe with :
- $\Phi = 20 \mu\text{m}$
  - Pitch =  $100\mu\text{m}$
  - ➔ **High aspect ratio 50 : 1**

Gold has been properly deposited over the whole surface and length of the holes




## Study of Heat Affected Zone (HAZ) : Crystallographic information (Electron Back-Scattering Diffraction SEM)


3 x 3 holes matrix drilled into 500 μm of GaAs  
Chemical etch with Br in methanol

- No recrystallisation and no amorphous area around laser drilled holes
- GaAs stays single crystal  $\langle 100 \rangle$

**No critical HAZ : HAZ suppressed by chemical treatment**




# 3D CdTe detector




Stabilization of the peak channel in both bias configuration above -10 V → full charge collection

Square-shaped bias configuration slightly more efficient than diamond-shaped bias configuration




# 3D CdTe detector



Stabilization of counts number  
for both bias configuration  
In the -20 to -50 V range

At -20V, ratio of total counts = 1.7



Demonstration of photon counting and  $^{241}\text{Am} / ^{57}\text{Co}$  peak energy discrimination using a non-optimized design  
(electrode diameter and pitch have to be optimized)

# 3D CdTe detector

