



# SiPM arrays and miniaturized readout electronics for compact imaging gamma camera

*N.Dinu<sup>1</sup>, T. Imando<sup>1</sup>, A. Nagai<sup>1</sup>, L. Pinot<sup>2</sup>, V. Puill<sup>1</sup>, S. Callier<sup>3</sup>, B. Janvier<sup>2</sup>,  
C. Esnault<sup>2</sup>, Y. Charon<sup>2</sup>, L. Raux<sup>3</sup>, V. Vandebussche<sup>2</sup>, M-A. Verdier<sup>2</sup>, L. Menard<sup>2</sup>*

<sup>1</sup>Laboratory of Linear Accelerator, IN2P3, CNRS, Orsay, France

<sup>2</sup>Laboratory of Imaging and Modelisation in Neurobiology and Cancerology, IN2P3, CNRS,  
Orsay, France


<sup>3</sup>Omega Microelectronics group, CNRS, Palaiseau, France

# OUTLINE

- Introduction
  - Technique of nuclear imaging
  - Radio-guided surgery
  - Per-operative detection systems
- Compact imaging gamma camera - MAGICS
  - Principle
  - SiPM arrays
  - Read-out electronics
  - Characteristics
- Conclusions


# Technique of nuclear imaging

## • Principle



### Marking

- radio-pharmaceutical product
- radioactive isotopes
  - $^{99m}\text{Tc}$ ,  $^{123}\text{I}$ ,  $^{201}\text{Tl}$ ,  $^{18}\text{F}$ ,  $^{11}\text{C}$
  - emitters  $\gamma$ ,  $\beta^+$





### Detection

- $\gamma$  camera, tomographs

## • Detection systems

### Cancer diagnostic



### Cancer therapy



# Cancer therapy

- Radio-guided surgery

- Surgical practice for different cancer pathologies

- protocol of sentinel lymphatic node mapping (breast, melanoma, prostate, uterus)
    - colorectal and thyroid tumors
    - neuro-endocrines tumors and their metastatic disseminations
    - subclinical tumors (non-palpable) of breast
    - bone lesions (osteoma)
    - parathyroid adenoma

- Per-operative detection systems

- Gamma probes

- ergonomic shape (pencil)
    - 1-2 cm diameter, 10-20 cm length
    - sound signal proportional to counting rate

- Gamma imaging cameras

- cover larger area: 10-100 cm<sup>2</sup>
    - give the spatial distribution of the radio-tracer
    - improve signal to noise ratio



POCI/TRECAM, France




Sentinella, Spain

# Compact imaging gamma camera

- High resolution hand-held radiation detector for therapeutic purposes



## MAGICS imaging camera




*Collaboration IMNC, LAL, Hôpital Lariboisière*

## Detection system requirements in surgical conditions



- reduced size and weight
- versatility of readout electronics
- high spatial resolution and sensitivity

## SI 1828-3344M Hamamatsu HPK

- 4x4 monolithic SiPM array
- mounted on a SMD package
- Each SiPM = one readout channel:
  - 3x3 mm<sup>2</sup>, 3600  $\mu$ cells, each  $\mu$ cell - 50x50  $\mu$ m<sup>2</sup>




# i-v measurements of monolithic SiPM arrays




Keithley 2611

Hi


Lo



I vs.  $V_{Bias}$




I vs.  $V_{Bias}$






# SiPM characteristics uniformity




$V_{BD}$  MAGICS camera @ 20°C




More details on SiPM measurements & characteristics:  
N. Dinu, HDR, LAL-13-192

$I_{dark} @ V_{BIAS} = 72.5V$




$I_{dark} @ \text{overvoltage} = 1V$



$$I_{post-BD} \sim q \cdot G \cdot DCR$$

Nicoleta Dinu, LAL, Orsay, France


# Elementary module of MAGICS camera



*T. Ait Imando et al, PoS 2012*

*Thanks to SERDI (LAL) for technical help*


Elementary module  
Field of view:  $\sim 8 \text{ cm}^2$



**Board 1:**  
4 (2x2) SiPM arrays  
64 channels

**Board 2:**  
2 EASIROC chips  
64 readout channels  
2 ADC 12 bits

**Board 3:**  
ALTERA cyclone III FPGA  
FTDI FT2232H (USB, 2.0 Hi-speed, 440MBit/s)  
DC/DC converter for SiPM bias






# Characteristics of read-out electronics

## EASIROC chip

- 32-channels fully analog front-end readout
- 8-bit DAC (0-2.5V) for individual SiPM gain adjustment
- energy measurement from 160 fC to 320 pC
  - 1 to 2000 pe @ SiPM gain of  $10^6$
  - variable gain pre-amplifier tuned to 4 bits
  - variable shaping time from 25 to 175 ns
  - 2 multiplexed analog outputs (high gain, low gain)
  - 1 pe signal/noise ratio  $\approx 9$
- Low power consumption
  - 4.84 mW/channel, 155 mW/chip


## Channels uniformity for HG output

$$\sigma < 1\%$$




## Channels linearity for HG output

$$\sigma < 1\% \text{ from } 160 \text{ fC to } 100 \text{ pC}$$




## SiPM gain correction vs temperature using input DAC



# MAGICS camera

4 elementary modules  
256 SiPM's = 256 readout channels

## Architecture




## Mechanics

- alignment and assembling


## Software

- boards driving, data acquisition and treatment

## 4 elementary modules



## MAGICS camera final view




**Dimensions:** 8.3 x 8.3 x 8.35 cm<sup>3</sup>

**Weight:** 1.2 kg

**Field of view :** 5.1x5.1 cm<sup>2</sup>

# Characteristics of MAGICS camera

Sources of  $^{241}\text{Am}$  &  $^{57}\text{Co}$



Energy resolution:

- 9.5% @ 122 keV

Experimental conditions:



- $\text{LaBr}_3(\text{Ce})$  – 6 mm thickness
- Sources of  $^{57}\text{Co}$  (122 keV) and  $^{241}\text{Am}$  (60 keV)
- central collimation hole of 4 mm diameter
- $V_{\text{BIAS}} = 75.52\text{V}$ ,  $T=40^\circ\text{C}$

Spatial resolution:

- 1.13 mm @ 122 keV

Experimental conditions:

- $\text{LaBr}_3(\text{Ce})$  – 6 mm thickness
- $^{57}\text{Co}$  (122 keV)
- 1 mm irradiation spots spaced 5 mm apart
- $V_{\text{BIAS}} = 75.52\text{V}$ ,  $T=40^\circ\text{C}$
- Levenberg-Marquard reconstruction algorithm



# Conclusions

- **MAGICS**

- Miniaturized hand-held gamma camera
- Therapeutic purposes: assisting the surgeon on locating and removing tumors

- **MAGICS main components**

- 256 SiPM's: 4x4 monolithic arrays, 4x4 SiPM's each array
- 8 EASIROC chips for analogue front-end read-out
- Dedicated miniaturized readout electronics boards for data acquisition

- **MAGICS main characteristics**

- Geometry:
  - $8.3 \times 8.3 \times 8.35 \text{ cm}^3$
  - *Field of view:  $5.1 \times 5.1 \text{ cm}^2$*
- Weight: 1.2 kg
- Detection:
  - Energy resolution of 9.5% @ 122 keV
  - Spatial resolution of less than 1 mm @ 122 keV on the whole field of view

- **Future improvements**

- Decrease the lateral dimensions and weight by using non-hygroscopic scintillator (GaGG)
- Improve the reconstruction algorithm






# Additional slides

# Comparison of MAGICS and TRECaM cameras

## MAGICS camera



Weight: 1.2 kg

## TRECaM camera based on MAPMT



Weight: 2.2 kg

MAGICS energy resolution: 9.5% @ 122 keV

TRECaM energy resolution: 12.9% @ 122 keV

MAGICS spatial resolution: less than 1 mm @ 122 keV

TRECaM spatial resolution: 1.36 mm @ 122 keV

# Uniformity adjustments

- Uniform illumination of SiPM's by green LED
- Standardization of responses by DACs adjustments

